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
Reference

1. A model for sidewall adhesions in the rabbit: Reduction by an absorbable barrier [Michael P. Diamond, Cary B. Linskv, Tim Cunningham, Barry Constantine, Gere S. Dizerega and Alan H. Decherney - MICROSURGERY 8:197-200 1987]
2. Foreign Material in Postoperative Adhesions. [Luijendijk and Others - ANN. SURG. Vol. 223, No. 3, March 1996]
3. Evaluation of the behaviour of three synthetic meshes for abdominal defects repair in rabbits [S.Cuccommarino, O.Capitani, M.Zanatta, S.Valentini, P.A.Crepaldi, R.Lamberti, J.J.Martinez Rodriguez - EHS 30 Congress Sevilla (Spain) Posters session VI. Miscelanea anatomy, biomaterials, basic research, socioeconomic asp, 7-10 May 2008]
4. Pathophysiologie und Klassifikation von Adhiisionen [Zuhlke HV, Lorenz EMP, Straub EM, Savva V. - Langenbecks Arch Chir Suppl II (Kongressbericht 1990): 1009-1016]
5. L'ernia addominale per via laparoscopica con protesi bimatereale Relimesh [G. Pozzo, G.M. Amerio, R. Bona, E. Castagna, C. Moro, U. Parisi, V. Sorisio, B.M.Marino, SOC Chirurgia Generale, Ospedale "Cardinal Massaia", ASTI, Italy - XX9° ACOI nazionale, Paestum (Italy) Chirurgia, tra mito e scienza, Communications 2, 26-29 May 2010]

Codici - Codes	Forma-Shape	Dimensioni - Dimensions	Pezzi/cf. - Units/box
PEP1114-O	Oval	cm 11 x 14	1
PEP1418-O	Oval	cm 14 x 18	1
PEP2025-O	Oval	cm 20 x 25	1
PEP2535-O	Oval	cm 25 x 35	1
PEP12-R	Round	Ø cm 12	1
PEP3030	Square	cm 30 x 30	1

A comparative experimentation has been carried out with animals, in order to evaluate the in vivo behaviour of Relimesh and its antagonistic products, with a particular regard to adhesion formation and inflammatory reactions.

Adhesions were classified according to Zuhlke⁴ scale and Relimesh showed a good behaviour in the intraperitoneal use, with a low adhesion formation.

■ Implantation of Relimesh in a rabbit by scraping the visceral peritoneum^{1,3}.

■ Explantation of the prosthesis after 48 hours: Relimesh shows no adhesions, so a classification according to Zuhlke is not possible.

■ Explantation of the prosthesis after 4 weeks: Relimesh shows a significant neoperitonealization and mild adhesions (grade II of Zuhlke scale) being equal or even lower comparable to those generated with the antagonist products evaluated in the same study. Moreover, the statistical study has demonstrated as Relimesh exhibits a stabler behavior, relatively to the interactions prosthesis - peritoneum, of the evaluated antagonist products^{2,3}.

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
- The microporosity of the ePTFE film is such as to minimize the risk of tissue adhesions and give a permeability which avoids the possible colonization of pathogenic bacteria. The softness of the coupled mesh guarantees a high comfort to the patient.

Relimesh® is

the prosthesis for the treatment of hernias by intraperitoneal way, which has been developed at Herniamesh laboratories and realized by coupling a light monofilament Polypropylene (PP) mesh with an expanded Polytetrafluoroethylene (ePTFE) film.


The ePTFE side has a very smooth, low-porosity surface, which is capable to minimize tissue adhesions, while the macroporous structure of the PP mesh stimulates a quick fibroblastic response, with a fast fixing to the host tissue as a result. Relimesh is patent pending product.

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Relimesh® is suitable for

hernias and abdominal wall defects such as incisional hernia, with open or laparoscopic approach. The easy identification of the two sides allows a correct positioning of the prosthesis, which is essential for a successful implantation: the PP side, with its rough surface, is designed for contact with the abdominal wall, while the smooth, white ePTFE side is specifically intended for visceral contact.

Lightness, resistance and reliability

Relimesh, with only 0.55 mm in thickness, provides values of tensile and burst strenght which are superior to the loads which could develop in human body under conditions of maximum tension, and gives a response which is compatible with normal physiological reactions as well as anatomical and biomechanical dynamics of abdominal wall.

Relimesh® can be

- cut to the required dimension and shape, without the prosthesis is affected by delamination or structural alterations. This thanks to the special mesh knitting and the exclusive process of combination between the two materials;
- easily rolled up for its introduction into the abdominal cavity by means of a trocar;
- laid easily on the anatomical site.